



# EYE-INSPIRED™ DESIGN ACUVUE® BRAND PARAMETER GUIDE.

# ACUVUE®

PACKAGE IMAGE	PRODUCT NAME	MATERIAL	WATER CONTENT	POWERS	DK/T	BASE CURVE/S	DIAMETER	CENTRE THICKNESS	VISIBILITY TINT	PACK SIZE	MODALITY	UV BLOCKING
	<b>ACUVUE® OASYS® with HYDRACLEAR® Plus</b>	Senofilcon A Silicone Hydrogel with Internal Wetting Agent	38%	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)	147	8.4, 8.8	14mm	0.07mm @ -3.00D	Light Blue	6 Pack	2 Week Daily Wear 1 Week Extended Wear	UV class I
	<b>ACUVUE® OASYS® for ASTIGMATISM</b>	Senofilcon A Silicone Hydrogel with Internal Wetting Agent	38%	+6.00D to -6.00D (0.25D steps) -6.50D to -9.00D (0.50D steps) Cyls: -0.75, -1.25, -1.75, -2.25, -2.75 Axes: 10° to 180° in 10° steps	129	8.6	14.5mm	0.08mm @ -3.00D	Light Blue	6 Pack	2 Week Daily Wear 1 Week Extended Wear	UV class I
	<b>ACUVUE OASYS® Brand Contact Lenses with HydraLuxe™</b>	Senofilcon A Silicone Hydrogel with Internal Wetting Agent	38%	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)	121	8.5, 9.0	14.3mm	0.085mm @ -3.00D	Light Blue	30 Pack & 90 Pack	Single Use Daily Wear	UV class I
	<b>ACUVUE® Oasys with Transitions™ Light Intelligent Technology™</b>	Senofilcon A (38%)	38%	-6.00D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps) -6.50D to -12.00D (0.50D steps) Not available in Plano	121	8.4	14mm	0.085mm @ -3.00D	Yes	6 Pack	2 Week Daily Wear	UV class I
	<b>1-DAY ACUVUE® MOIST Brand MULTIFOCAL Contact Lenses</b>	Etafilcon A hydrogel lens with Internal Wetting Agent	58%	Spherical Power -9.00D to +6.00D (0.25D steps) Adds: +0.75D to +1.25D (LOW) +1.50D to +1.75D (MID) +2.00D to +2.50D (HGH)	25.5	8.4	14.3mm	0.084mm @ -3.00D (Varies with power)	Light Blue	30 Pack	Single Use Daily Wear	UV class II
	<b>1-DAY ACUVUE® MOIST®</b>	Etafilcon A hydrogel lens with Internal Wetting Agent	58%	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps)	25.5	8.5, 9.0	14.2mm	0.084mm @ -3.00D	Light Blue	30 Pack & 90 Pack	Single Use Daily Wear	UV class II
	<b>1-DAY ACUVUE® MOIST® for ASTIGMATISM</b>	Etafilcon A hydrogel lens with Internal Wetting Agent	58%	Plano to -6.00D (0.25D steps) Cyls: -0.75, -1.25, -1.75 Axes: 10° to 180° in 10° steps Cyls: -2.25 Axes: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180° -6.50D to -9.00D (0.50D steps), +0.25D to +4.00D (0.25D steps) Cyls: -0.75, -1.25, -1.75 Axes: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°	25.5	8.5	14.5mm	0.09mm @ -3.00D	Light Blue	30 Pack	Single Use Daily Wear	UV class II
	<b>ACUVUE® 2</b>	Etafilcon A	58%	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)	25.5	8.3, 8.7	14mm	0.084mm @ -3.00D	Light Blue	6 Pack	2 Week Daily Wear	UV class II

All ACUVUE® Brand Contact Lenses have class 1 or Class 2 UV Blocking to help provide protection against transmission of harmful UV radiation to the cornea and into the eye. UV absorbing contact lenses are NOT substitutes for protective UV absorbing eyewear such as UV absorbing goggles or sunglasses because they do not completely cover the eye and surrounding area. You should continue to use UV absorbing eyewear as directed by your Eye Care Professional. ACUVUE®2®, 1-DAY ACUVUE® MOIST, ACUVUE OASYS®, HYDRACLEAR® PLUS, 1-DAY ACUVUE® TruEye®, LACREON® and EYE-INSPIRED™ DESIGN are trademarks of Johnson & Johnson Medical (Pty) Ltd. © Johnson & Johnson Medical (Pty) Ltd. 2019. Transitions™ is a registered trademark, the Transitions™ logo, Transitions™ Light Intelligent Technology™ and Transitions Light Intelligent™ Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. ©Johnson & Johnson Medical (Pty) Ltd. 2019. MOS126\_19